



I
N
F
O
R
M
A
T
I
O
N

Bubble Pressure Tensiometer BP100



Printing – bonding – spraying – cleaning: processes such as these require surfactants that are already effective in fractions of a second – and an instrument that can record their effectiveness. The Bubble Pressure Tensiometer BP100 determines the dynamic surface tension as a function of the surface age. Thus the instrument provides information about wetting and drop formation in rapid processes.

- **Measurement of the dynamic surface tension at constant or variable surface age**
- **Wide range of surface age:
5 ms (high dynamic) up to 100 s (almost static)**
- **Fully automatic measuring process by software controlled immersion of the capillary**
- **Determination of the adsorption and diffusion coefficient**
- **Integrated compressor -
no compressed air supply necessary**



Bubble Pressure Measurement with the BP100

A software-controlled flow of air emerging from a capillary produces air bubbles in the sample. A high-precision pressure sensor determines the maximum pressure during bubble formation, from which the surface tension is calculated. Thanks to its built-in compressor the instrument does not need an external compressed air connection. Thermostatted measurements are possible without any problems.

Scans with an almost unlimited resolution over a wide range of surface ages are controlled automatically by the LabDesk software. In addition to the time-dependent surface tension, the measurement also supplies the equilibrium value according to Hua & Rosen. Diffusion and adsorption coefficients can be calculated from surfactant concentration series – these are important parameters where the mobility of surfactants is concerned. Measurements at constant bubble ages allow further insights, e.g. for concentration or temperature comparisons.

After a minimal preparation period the measurement is carried out fully automatically up to data output in the form of a plot. With a mouse-click measuring parameters, results and evaluations appear in a comprehensive report. Thanks to the LabDesk™ software platform for all KRÜSS tensiometers data from other instruments can be included.



Applications

- Surface-active agent development
- Optimization of spraying processes
- Development of detergents and cleansing agents
- Optimization of coating and printing processes
- Electroplating bath concentration control

Technische Daten

Measuring range:

- SFT ⁽¹⁾: 10 to 100 mN/m
- Temperature: -10 to 100 °C

Resolution:

- SFT: ± 0.01 mN/m
- Temperature: ± 0.01°

Surface age ⁽¹⁾: 5 ms - 100 s

Interfaces: USB, RS232C

Power consumption: max. 30 W

Power supply: 100 – 240 V AC; 47 – 63 Hz

Weight: 15 kg

Dimensions: 300 x 540 x 320 mm (WxHxL)

(1): depending on capillary and liquid

Measuring Methods

- Measurement of surface tension at constant surface age
- Measurement of surface tension at any surface ages from 5 to 100,000 ms

Accessories

- Various thermostats
- Capillaries made of glass or PTFE
- Coating Kit to renew the hydrophobicity of the glass capillaries
- Thermostat vessel with built-in magnetic stirrer

Technical specifications are subject to change without notice.



<http://www.kruss.de>

KRÜSS GmbH
Wissenschaftliche Laborgeräte
Borsteler Chaussee 85-99a
D-22453 Hamburg
Tel: +49 - 40 - 51 44 01 - 0
Fax: +49 - 40 - 51 44 01 - 98
E-Mail: info@kruss.de

KRÜSS GmbH
38/40 Avenue Jean Jaurès
F-91120 Palaiseau
Tel: +33 - 1 - 60 14 94 94
Fax: +33 - 1 - 60 14 95 48
E-Mail: info@kruss.fr

KRÜSS USA
1020 Crews Road, Suite K
Matthews, NC 28105

Tel: +1 - 704 - 847 8933
Fax: +1 - 704 - 847 9416
E-Mail: info@kruss-usa.com